



Polishine®

PLL

Designed for two steps polishing of Lithium Disilicate (e.g. IPS e.max®) restorations.

Lithium Disilicate System Including Autoclavable Bur Block

KIT PLL - Lithium Disilicate System

2 Steps Set

1st - Silk Finish: Blue

2nd - Hi Gloss Polishing: Yellow

Zirconium Oxide (ZrO₂) ceramic material for finishing and high gloss polish of Lithium Disilicate

Pack of: 3+3

Autoclavable 135°C

RA SS shanks

1st step rec. RPM: 7,000-12,000. Max. 20,000


2nd step rec. RPM: 7,000-10,000. Max. 20,000



Order No. kit PLL

+FREE BUR BLOCK

Lithium Disilicate Polishers



	1st - Silk Finish			2nd - Hi-Gloss Polishing		
	disc	cone	cup	disc	cone	cup
Item #	PLL 101SF	PLL 510SF	PLL 675SF	PLL 101HG	PLL 510HG	PLL 675HG
Dimensions Ø X L1	10 x 0.7	5 x 10	6.5 x 9.5	10 x 0.7	5 x 10	6.5 x 9.5
Color	Turquoise	Turquoise	Turquoise	Yellow	Yellow	Yellow
Grit	■ M	■ M	■ M	■ F	■ F	■ F
Recommended RPM	7,000 - 12,000	7,000 - 12,000	7,000 - 12,000	7,000-10,000	7,000-10,000	7,000-10,000

* IPS e.max is a registered trademark of IVOCCLAR VIVADENT, Inc.